

DR 4K Camera



The DR MEDICAL camera system is one high-resolution camera for capturing still and video images in endoscopic and general surgical applications.


The system also contains an internal light source, via a turret optical fiber adapter has, on which the different light guides fit.

This DR 4K camera head is autoclavable.

Advantages

- 4K UHD resolution
- Zoom, freeze, fog reduction
- Live videos and still images
- Archiving and documenting pictures and videos

Color space



Technical parameters

DR 4K camera with image processor, LED light source & data acquisition

Ref. No. 85-10-4010

Voltage range	100-240 V 50-60 Hz 400 VA
Video outputs	DisplayPort (4K): 3840x2160, progressive scanning DVI (1080P): 1920x1080, progressive scanning, SMPTE 425 1920x1080, progressive scanning
Connections	DVI Out, Displayport, DVI In, Still Out, Rec Out, Foot In, Mic In, Tablet, USB 2.0, peripherals, Network, USB 3.0, AC input, equipotential bonding ground connection
Monitor	Touchscreen
Fiber optic adapter	Revolver for ACMI, Storz, Wolf, Olympus
White balance range	3000 to 7500 K
Dimensions (W x H x D)	305 x 165 x 368 mm
Weight	5.9 kg



DR 4K camera head

Ref. No. 85-10-4100

Casing	With 4 programmable pushbuttons
Camera head	Autoclavable and for STERRAD® cleanings
Diameter	40 mm
Length	89 mm
Weight	510 g



DR Tablet

Ref. No. 85-10-1000

Cleaning	Wipe disinfection with commercially available agents
Color	Black
Dimensions (W x H x D)	254 x 155 x 8 mm
Weight	530 g



Legal Manufacture: Santa Barbara Imaging Systems / 340 Storke Rd, Suite 101 / Goleta, CA 93117 USA

Version001_EN

DR MEDICAL AG
Ziegelmatteareal, Bau 6
CH- 4503 Solothurn
SWITZERLAND

T. + 41 32 621 60 01
F. + 41 32 621 60 02
info@dr-medical.ch
www.dr-medical.ch

SWITZERLAND
DR MEDICAL 
for your arthroscopy solution