

DR Bipolar Arthroscopy Electrode System



The DR Bipolar Arthroscopy Electrode System is equipped with the most modern technology and software.

The electrodes serve to cut, vaporize, coagulate, and thermally alter soft tissue.

The generator was specially developed for arthroscopy.

Bipolar arthroscopy system with the newest technology

Advantages


- High-tech bipolar mode for work in sodium chloride solution
- Universally deployable across the bipolar and monopolar spectrum of arthroscopy
- Automatic control of parameters / performance during activation
- Instrument detection with default settings that preclude excessively high values

Uses



Shoulder joints

Acromioplasty
Partial synovectomy
Cartilage smoothing
Tenotomy



Knee joints

Partial meniscal resection
Plica synovialis resection
Partial synovectomy
Lateral release



Hip joints

Synovectomy
Cartilage smoothing
Acromioplasty
Partial labral resection

DR High-End RF Products

DR Bipolar Arthroscopy Electrode Generator

Ref. No. 897-6

Display	Thin-film-transistor 7" wide video graphics array
Technology	Touchscreen
Spectrum	Automatic performance control
Dimensions (w x h x d)	305 x 141 x 385 mm
Weight	7 kg



DR Handpiece with Connecting Cable

Ref. No. 92888-6

Tube connector for aspiration
Four-metre connecting cable
Guaranteed for three hundred procedures



DR Handle with Two Buttons

Ref. No. 92889-6

DR handle for 4-mm-shaft electrodes with integrated cable and two buttons.
Guaranteed for three hundred procedures



DR Generator Foot Pedal

Ref. No. 897-7

Foot pedal (two pedals) for the DR Bipolar Arthroscopy Electrode Generator.



DR Vaporisation Electrode with Aspirator, 90°, 150 mm long

Ref. No. 890-1s

Sterile packaging, box of ten



DR Vaporisation Electrode with Aspirator, 45°, 150 mm long

Ref. No. 890-6s

Sterile packaging, box of ten



DR Vaporisation Electrode with Aspirator, 45°, 150 mm long

Ref. No. 890-5s

Sterile packaging, box of ten



DR Vaporisation Electrode with Aspirator, 45°, 180 mm long

Ref. No. 891-2s

Sterile packaging, box of ten



DR Vaporisation Electrode with Aspirator, flexible to 40°, 150 mm long

Ref. No. 890-7s

Sterile packaging, box of ten



* Further electrode varieties are available; please enquire.

Version001_EN