

DR Shaver Blades and Burs for the ARTHREX Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	023.801.230	DRAR
DR WAVE CUTTER	Ø 3.0	85	023.301.231	DRAR
DR WAVE CUTTER	Ø 3.0	130	023.301.230	DRAR
DR WAVE CUTTER	Ø 4.0	130	023.101.230	DRAR
DR WAVE CUTTER	Ø 5.0	130	023.201.230	DRAR
DR PIT BULL	Ø 3.0	130	023.301.240	DRAR
DR PIT BULL	Ø 4.0	130	023.101.240	DRAR
DR PIT BULL	Ø 5.0	130	023.201.240	DRAR
DR PIT BULL	Ø 7.0	130	023.701.240	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.236	DRAR
DR BENDING PIT BULL	Ø 3.0	130	023.301.340	DRAR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.241	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.242	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.250	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.251	DRAR
DR BARREL BUR	Ø 4.0	130	133.100.220	DRAR
DR BARREL BUR	Ø 5.5	130	133.200.220	DRAR
DR ROUND BUR	Ø 4.0	130	123.100.220	DRAR
DR ROUND BUR	Ø 5.5	130	123.200.220	DRAR
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	023.801.030	DRAR
DR WAVE CUTTER	Ø 3.0	85	023.301.030	DRAR
DR WAVE CUTTER	Ø 3.0	130	023.301.031	DRAR
DR WAVE CUTTER	Ø 4.0	130	023.101.030	DRAR
DR WAVE CUTTER	Ø 5.0	130	023.201.030	DRAR
DR PIT BULL	Ø 3.0	85	023.301.047	DRAR
DR PIT BULL	Ø 3.0	130	023.301.040	DRAR
DR PIT BULL	Ø 4.0	130	023.101.040	DRAR
DR PIT BULL	Ø 5.0	130	023.201.040	DRAR
DR PIT BULL	Ø 7.0	130	023.701.040	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.036	DRAR
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	023.201.043	DRAR
DR PIT BULL	Ø 4.0	155	023.101.037	DRAR
DR PIT BULL	Ø 5.0	155	023.201.037	DRAR
DR BENDING PIT BULL	Ø 3.0	130	023.301.140	DRAR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.041	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	023.101.042	DRAR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.050	DRAR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	023.101.051	DRAR
DR BARREL BUR	Ø 4.0	130	133.100.020	DRAR
DR BARREL BUR	Ø 5.5	130	133.200.020	DRAR
DR BARREL BUR	Ø 4.0	174	133.100.021	DRAR
DR ROUND BUR	Ø 4.0	130	123.100.020	DRAR
DR ROUND BUR	Ø 5.5	130	123.200.020	DRAR
DR ROUND BUR	Ø 4.0	174	123.100.021	DRAR
DR ROUND BUR	Ø 5.5	174	123.200.021	DRAR
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	023.101.528	DRAR
DR WAVE CUTTER	Ø 5.0	130	023.201.528	DRAR
DR PIT BULL	Ø 4.0	130	023.101.538	DRAR
DR PIT BULL	Ø 5.0	130	023.201.538	DRAR
DR PIT BULL	Ø 7.0	130	023.701.538	DRAR
DR PIT BULL	Ø 4.0	155	023.101.535	DRAR
DR PIT BULL	Ø 5.0	155	023.201.535	DRAR
DR SYNOVIA PIT BULL	Ø 4.0	130	023.101.534	DRAR
DR BARREL BUR	Ø 4.0	130	133.100.519	DRAR
DR BARREL BUR	Ø 5.5	130	133.200.519	DRAR
DR ROUND BUR	Ø 4.0	130	123.100.519	DRAR
DR ROUND BUR	Ø 5.5	130	123.200.519	DRAR

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)



DR PIT BULL®

(three double blades, large window)




DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRAR DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by ARTHREX. The firm Arthrex Corporation Inc. and DR MEDICAL AG are independent companies.

Version001_EN