

DR Shaver Blades and Burs for the DR MEDICAL Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	029.801.230	DRVI
DR WAVE CUTTER	Ø 3.0	85	029.301.231	DRVI
DR WAVE CUTTER	Ø 3.0	130	029.301.230	DRVI
DR WAVE CUTTER	Ø 4.0	130	029.101.230	DRVI
DR WAVE CUTTER	Ø 5.0	130	029.201.230	DRVI
DR PIT BULL	Ø 3.0	130	029.301.240	DRVI
DR PIT BULL	Ø 4.0	130	029.101.240	DRVI
DR PIT BULL	Ø 5.0	130	029.201.240	DRVI
DR PIT BULL	Ø 7.0	130	029.701.240	DRVI
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	029.101.236	DRVI
DR BENDING PIT BULL	Ø 3.0	130	029.301.340	DRVI
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	029.101.241	DRVI
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	029.101.242	DRVI
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	029.101.250	DRVI
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	029.101.251	DRVI
DR BARREL BUR	Ø 4.0	130	139.100.220	DRVI
DR BARREL BUR	Ø 5.5	130	139.200.220	DRVI
DR ROUND BUR	Ø 4.0	130	129.100.220	DRVI
DR ROUND BUR	Ø 5.5	130	129.200.220	DRVI
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	029.801.030	DRVI
DR WAVE CUTTER	Ø 3.0	85	029.301.030	DRVI
DR WAVE CUTTER	Ø 3.0	130	029.301.031	DRVI
DR WAVE CUTTER	Ø 4.0	130	029.101.030	DRVI
DR WAVE CUTTER	Ø 5.0	130	029.201.030	DRVI
DR PIT BULL	Ø 3.0	85	029.301.047	DRVI
DR PIT BULL	Ø 3.0	130	029.301.040	DRVI
DR PIT BULL	Ø 4.0	130	029.101.040	DRVI
DR PIT BULL	Ø 5.0	130	029.201.040	DRVI
DR PIT BULL	Ø 7.0	130	029.701.040	DRVI
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	029.101.036	DRVI
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	029.201.043	DRVI
DR PIT BULL	Ø 4.0	155	029.101.037	DRVI
DR PIT BULL	Ø 5.0	155	029.201.037	DRVI
DR BENDING PIT BULL	Ø 3.0	130	029.301.140	DRVI
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	029.101.041	DRVI
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	029.101.042	DRVI
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	029.101.050	DRVI
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	029.101.051	DRVI
DR BARREL BUR	Ø 4.0	130	139.100.020	DRVI
DR BARREL BUR	Ø 5.5	130	139.200.020	DRVI
DR BARREL BUR	Ø 4.0	174	139.100.021	DRVI
DR ROUND BUR	Ø 4.0	130	129.100.020	DRVI
DR ROUND BUR	Ø 5.5	130	129.200.020	DRVI
DR ROUND BUR	Ø 4.0	174	129.100.021	DRVI
DR ROUND BUR	Ø 5.5	174	129.200.021	DRVI
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	029.101.528	DRVI
DR WAVE CUTTER	Ø 5.0	130	029.201.528	DRVI
DR PIT BULL	Ø 4.0	130	029.101.538	DRVI
DR PIT BULL	Ø 5.0	130	029.201.538	DRVI
DR PIT BULL	Ø 7.0	130	029.701.538	DRVI
DR PIT BULL	Ø 4.0	155	029.101.535	DRVI
DR PIT BULL	Ø 5.0	155	029.201.535	DRVI
DR SYNOVIA PIT BULL	Ø 4.0	130	029.101.534	DRVI
DR BARREL BUR	Ø 4.0	130	139.100.519	DRVI
DR BARREL BUR	Ø 5.5	130	139.200.519	DRVI
DR ROUND BUR	Ø 4.0	130	129.100.519	DRVI
DR ROUND BUR	Ø 5.5	130	129.200.519	DRVI

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)




DR PIT BULL®

(three double blades, large window)




DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®