

DR Shaver Blades and Burs for the STORZ Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	026.801.230	DRSZ
DR WAVE CUTTER	Ø 3.0	85	026.301.231	DRSZ
DR WAVE CUTTER	Ø 3.0	130	026.301.230	DRSZ
DR WAVE CUTTER	Ø 4.0	130	026.101.230	DRSZ
DR WAVE CUTTER	Ø 5.0	130	026.201.230	DRSZ
DR PIT BULL	Ø 3.0	130	026.301.240	DRSZ
DR PIT BULL	Ø 4.0	130	026.101.240	DRSZ
DR PIT BULL	Ø 5.0	130	026.201.240	DRSZ
DR PIT BULL	Ø 7.0	130	026.701.240	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.236	DRSZ
DR BENDING PIT BULL	Ø 3.0	130	026.301.340	DRSZ
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.241	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.242	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.250	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.251	DRSZ
DR BARREL BUR	Ø 4.0	130	136.100.220	DRSZ
DR BARREL BUR	Ø 5.5	130	136.200.220	DRSZ
DR ROUND BUR	Ø 4.0	130	126.100.220	DRSZ
DR ROUND BUR	Ø 5.5	130	126.200.220	DRSZ
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	026.801.030	DRSZ
DR WAVE CUTTER	Ø 3.0	85	026.301.030	DRSZ
DR WAVE CUTTER	Ø 3.0	130	026.301.031	DRSZ
DR WAVE CUTTER	Ø 4.0	130	026.101.030	DRSZ
DR WAVE CUTTER	Ø 5.0	130	026.201.030	DRSZ
DR PIT BULL	Ø 3.0	85	026.301.047	DRSZ
DR PIT BULL	Ø 3.0	130	026.301.040	DRSZ
DR PIT BULL	Ø 4.0	130	026.101.040	DRSZ
DR PIT BULL	Ø 5.0	130	026.201.040	DRSZ
DR PIT BULL	Ø 7.0	130	026.701.040	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.036	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	026.201.043	DRSZ
DR PIT BULL	Ø 4.0	155	026.101.037	DRSZ
DR PIT BULL	Ø 5.0	155	026.201.037	DRSZ
DR BENDING PIT BULL	Ø 3.0	130	026.301.140	DRSZ
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.041	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	026.101.042	DRSZ
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.050	DRSZ
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	026.101.051	DRSZ
DR BARREL BUR	Ø 4.0	130	136.100.020	DRSZ
DR BARREL BUR	Ø 5.5	130	136.200.020	DRSZ
DR BARREL BUR	Ø 4.0	174	136.100.021	DRSZ
DR ROUND BUR	Ø 4.0	130	126.100.020	DRSZ
DR ROUND BUR	Ø 5.5	130	126.200.020	DRSZ
DR ROUND BUR	Ø 4.0	174	126.100.021	DRSZ
DR ROUND BUR	Ø 5.5	174	126.200.021	DRSZ
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	026.101.528	DRSZ
DR WAVE CUTTER	Ø 5.0	130	026.201.528	DRSZ
DR PIT BULL	Ø 4.0	130	026.101.538	DRSZ
DR PIT BULL	Ø 5.0	130	026.201.538	DRSZ
DR PIT BULL	Ø 7.0	130	026.701.538	DRSZ
DR PIT BULL	Ø 4.0	155	026.101.535	DRSZ
DR PIT BULL	Ø 5.0	155	026.201.535	DRSZ
DR SYNOVIA PIT BULL	Ø 4.0	130	026.101.534	DRSZ
DR BARREL BUR	Ø 4.0	130	136.100.519	DRSZ
DR BARREL BUR	Ø 5.5	130	136.200.519	DRSZ
DR ROUND BUR	Ø 4.0	130	126.100.519	DRSZ
DR ROUND BUR	Ø 5.5	130	126.200.519	DRSZ

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)




DR PIT BULL®

(three double blades, large window)



DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRSZ DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by STORZ.

The firm Karl Storz SE & Co KG and DR MEDICAL AG are independent companies.

Version001_EN