

DR Shaver Blades and Burs for the STRYKER Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	022.801.230	DRST
DR WAVE CUTTER	Ø 3.0	85	022.301.231	DRST
DR WAVE CUTTER	Ø 3.0	130	022.301.230	DRST
DR WAVE CUTTER	Ø 4.0	130	022.101.230	DRST
DR WAVE CUTTER	Ø 5.0	130	022.201.230	DRST
DR PIT BULL	Ø 3.0	130	022.301.240	DRST
DR PIT BULL	Ø 4.0	130	022.101.240	DRST
DR PIT BULL	Ø 5.0	130	022.201.240	DRST
DR PIT BULL	Ø 7.0	130	022.701.240	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.236	DRST
DR BENDING PIT BULL	Ø 3.0	130	022.301.340	DRST
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.241	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.242	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.250	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.251	DRST
DR BARREL BUR	Ø 4.0	130	132.100.220	DRST
DR BARREL BUR	Ø 5.5	130	132.200.220	DRST
DR ROUND BUR	Ø 4.0	130	122.100.220	DRST
DR ROUND BUR	Ø 5.5	130	122.200.220	DRST
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	022.801.030	DRST
DR WAVE CUTTER	Ø 3.0	85	022.301.030	DRST
DR WAVE CUTTER	Ø 3.0	130	022.301.031	DRST
DR WAVE CUTTER	Ø 4.0	130	022.101.030	DRST
DR WAVE CUTTER	Ø 5.0	130	022.201.030	DRST
DR PIT BULL	Ø 3.0	85	022.301.047	DRST
DR PIT BULL	Ø 3.0	130	022.301.040	DRST
DR PIT BULL	Ø 4.0	130	022.101.040	DRST
DR PIT BULL	Ø 5.0	130	022.201.040	DRST
DR PIT BULL	Ø 7.0	130	022.701.040	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.036	DRST
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	022.201.043	DRST
DR PIT BULL	Ø 4.0	155	022.101.037	DRST
DR PIT BULL	Ø 5.0	155	022.201.037	DRST
DR BENDING PIT BULL	Ø 3.0	130	022.301.140	DRST
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.041	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	022.101.042	DRST
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.050	DRST
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	022.101.051	DRST
DR BARREL BUR	Ø 4.0	130	132.100.020	DRST
DR BARREL BUR	Ø 5.5	130	132.200.020	DRST
DR BARREL BUR	Ø 4.0	174	132.100.021	DRST
DR ROUND BUR	Ø 4.0	130	122.100.020	DRST
DR ROUND BUR	Ø 5.5	130	122.200.020	DRST
DR ROUND BUR	Ø 4.0	174	122.100.021	DRST
DR ROUND BUR	Ø 5.5	174	122.200.021	DRST
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	022.101.528	DRST
DR WAVE CUTTER	Ø 5.0	130	022.201.528	DRST
DR PIT BULL	Ø 4.0	130	022.101.538	DRST
DR PIT BULL	Ø 5.0	130	022.201.538	DRST
DR PIT BULL	Ø 7.0	130	022.701.538	DRST
DR PIT BULL	Ø 4.0	155	022.101.535	DRST
DR PIT BULL	Ø 5.0	155	022.201.535	DRST
DR SYNOVIA PIT BULL	Ø 4.0	130	022.101.534	DRST
DR BARREL BUR	Ø 4.0	130	132.100.519	DRST
DR BARREL BUR	Ø 5.5	130	132.200.519	DRST
DR ROUND BUR	Ø 4.0	130	122.100.519	DRST
DR ROUND BUR	Ø 5.5	130	122.200.519	DRST

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)




DR PIT BULL®

(three double blades, large window)



DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRST DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by STRYKER. The firm Stryker Corporation Inc. and DR MEDICAL AG are independent companies.

Version001_EN