

DR Shaver Blades and Burs for the WOLF Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	025.801.230	DRWO
DR WAVE CUTTER	Ø 3.0	85	025.301.231	DRWO
DR WAVE CUTTER	Ø 3.0	130	025.301.230	DRWO
DR WAVE CUTTER	Ø 4.0	130	025.101.230	DRWO
DR WAVE CUTTER	Ø 5.0	130	025.201.230	DRWO
DR PIT BULL	Ø 3.0	130	025.301.240	DRWO
DR PIT BULL	Ø 4.0	130	025.101.240	DRWO
DR PIT BULL	Ø 5.0	130	025.201.240	DRWO
DR PIT BULL	Ø 7.0	130	025.701.240	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.236	DRWO
DR BENDING PIT BULL	Ø 3.0	130	025.301.340	DRWO
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.241	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.242	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.250	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.251	DRWO
DR BARREL BUR	Ø 4.0	130	135.100.220	DRWO
DR BARREL BUR	Ø 5.5	130	135.200.220	DRWO
DR ROUND BUR	Ø 4.0	130	125.100.220	DRWO
DR ROUND BUR	Ø 5.5	130	125.200.220	DRWO
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	025.801.030	DRWO
DR WAVE CUTTER	Ø 3.0	85	025.301.030	DRWO
DR WAVE CUTTER	Ø 3.0	130	025.301.031	DRWO
DR WAVE CUTTER	Ø 4.0	130	025.101.030	DRWO
DR WAVE CUTTER	Ø 5.0	130	025.201.030	DRWO
DR PIT BULL	Ø 3.0	85	025.301.047	DRWO
DR PIT BULL	Ø 3.0	130	025.301.040	DRWO
DR PIT BULL	Ø 4.0	130	025.101.040	DRWO
DR PIT BULL	Ø 5.0	130	025.201.040	DRWO
DR PIT BULL	Ø 7.0	130	025.701.040	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.036	DRWO
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	025.201.043	DRWO
DR PIT BULL	Ø 4.0	155	025.101.037	DRWO
DR PIT BULL	Ø 5.0	155	025.201.037	DRWO
DR BENDING PIT BULL	Ø 3.0	130	025.301.140	DRWO
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.041	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	025.101.042	DRWO
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.050	DRWO
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	025.101.051	DRWO
DR BARREL BUR	Ø 4.0	130	135.100.020	DRWO
DR BARREL BUR	Ø 5.5	130	135.200.020	DRWO
DR BARREL BUR	Ø 4.0	174	135.100.021	DRWO
DR ROUND BUR	Ø 4.0	130	125.100.020	DRWO
DR ROUND BUR	Ø 5.5	130	125.200.020	DRWO
DR ROUND BUR	Ø 4.0	174	125.100.021	DRWO
DR ROUND BUR	Ø 5.5	174	125.200.021	DRWO
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	025.101.528	DRWO
DR WAVE CUTTER	Ø 5.0	130	025.201.528	DRWO
DR PIT BULL	Ø 4.0	130	025.101.538	DRWO
DR PIT BULL	Ø 5.0	130	025.201.538	DRWO
DR PIT BULL	Ø 7.0	130	025.701.538	DRWO
DR PIT BULL	Ø 4.0	155	025.101.535	DRWO
DR PIT BULL	Ø 5.0	155	025.201.535	DRWO
DR SYNOVIA PIT BULL	Ø 4.0	130	025.101.534	DRWO
DR BARREL BUR	Ø 4.0	130	135.100.519	DRWO
DR BARREL BUR	Ø 5.5	130	135.200.519	DRWO
DR ROUND BUR	Ø 4.0	130	125.100.519	DRWO
DR ROUND BUR	Ø 5.5	130	125.200.519	DRWO

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)




DR PIT BULL®

(three double blades, large window)



DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRWO DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by WOLF. The firm Richard Wolf GmbH and DR MEDICAL AG are independent companies.

Version001_EN