

DR Shaver Blades and Burs for the Johnson & Johnson FMS DUO Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

| Description | Ø (mm) | L (mm) | Ref.No. | Hub |
|------------------------------------|--|--------|-------------|-----|
| Plastic blue hub | Single use instruments into sterile packaging with 5 pieces inside | | | |
| DR WAVE CUTTER | Ø 2.0 | 85 | 020.801.230 | FMS |
| DR WAVE CUTTER | Ø 3.0 | 85 | 020.301.231 | FMS |
| DR WAVE CUTTER | Ø 3.0 | 130 | 020.301.230 | FMS |
| DR WAVE CUTTER | Ø 4.0 | 130 | 020.101.230 | FMS |
| DR WAVE CUTTER | Ø 5.0 | 130 | 020.201.230 | FMS |
| DR PIT BULL | Ø 3.0 | 130 | 020.301.240 | FMS |
| DR PIT BULL | Ø 4.0 | 130 | 020.101.240 | FMS |
| DR PIT BULL | Ø 5.0 | 130 | 020.201.240 | FMS |
| DR PIT BULL | Ø 7.0 | 130 | 020.701.240 | FMS |
| DR PIT BULL, CONCAVE 15° | Ø 4.0 | 130 | 020.101.236 | FMS |
| DR BENDING PIT BULL | Ø 3.0 | 130 | 020.301.340 | FMS |
| DR BENDING TOOL | | | 020.101.103 | |
| DR SYNOVIA PIT BULL | Ø 4.0 | 130 | 020.101.241 | FMS |
| DR SYNOVIA PIT BULL, CONCAVE 15° | Ø 4.0 | 130 | 020.101.242 | FMS |
| DR SYNOVIA PIT BULL, CONCAVE 15° | Ø 4.0 | 174 | 020.101.250 | FMS |
| DR PIT BULL, CONCAVE 15° | Ø 4.0 | 174 | 020.101.251 | FMS |
| DR BARREL BUR | Ø 4.0 | 130 | 130.100.220 | FMS |
| DR BARREL BUR | Ø 5.5 | 130 | 130.200.220 | FMS |
| DR ROUND BUR | Ø 4.0 | 130 | 120.100.220 | FMS |
| DR ROUND BUR | Ø 5.5 | 130 | 120.200.220 | FMS |
| Metal hub silver | Single use instruments into sterile packaging with 5 pieces inside | | | |
| DR WAVE CUTTER | Ø 2.0 | 85 | 020.801.030 | FMS |
| DR WAVE CUTTER | Ø 3.0 | 85 | 020.301.030 | FMS |
| DR WAVE CUTTER | Ø 3.0 | 130 | 020.301.031 | FMS |
| DR WAVE CUTTER | Ø 4.0 | 130 | 020.101.030 | FMS |
| DR WAVE CUTTER | Ø 5.0 | 130 | 020.201.030 | FMS |
| DR PIT BULL | Ø 3.0 | 85 | 020.301.047 | FMS |
| DR PIT BULL | Ø 3.0 | 130 | 020.301.040 | FMS |
| DR PIT BULL | Ø 4.0 | 130 | 020.101.040 | FMS |
| DR PIT BULL | Ø 5.0 | 130 | 020.201.040 | FMS |
| DR PIT BULL | Ø 7.0 | 130 | 020.701.040 | FMS |
| DR PIT BULL, CONCAVE 15° | Ø 4.0 | 130 | 020.101.036 | FMS |
| DR PIT BULL, CONCAVE 15° | Ø 5.0 | 155 | 020.201.043 | FMS |
| DR PIT BULL | Ø 4.0 | 155 | 020.101.037 | FMS |
| DR PIT BULL | Ø 5.0 | 155 | 020.201.037 | FMS |
| DR BENDING PIT BULL | Ø 3.0 | 130 | 020.301.140 | FMS |
| DR BENDING TOOL | | | 020.101.103 | |
| DR SYNOVIA PIT BULL | Ø 4.0 | 130 | 020.101.041 | FMS |
| DR SYNOVIA PIT BULL, CONCAVE 15° | Ø 4.0 | 130 | 020.101.042 | FMS |
| DR SYNOVIA PIT BULL, CONCAVE 15° | Ø 4.0 | 174 | 020.101.050 | FMS |
| DR PIT BULL, CONCAVE 15° | Ø 4.0 | 174 | 020.101.051 | FMS |
| DR BARREL BUR | Ø 4.0 | 130 | 130.100.020 | FMS |
| DR BARREL BUR | Ø 5.5 | 130 | 130.200.020 | FMS |
| DR BARREL BUR | Ø 4.0 | 174 | 130.100.021 | FMS |
| DR ROUND BUR | Ø 4.0 | 130 | 120.100.020 | FMS |
| DR ROUND BUR | Ø 5.5 | 130 | 120.200.020 | FMS |
| DR ROUND BUR | Ø 4.0 | 174 | 120.100.021 | FMS |
| DR ROUND BUR | Ø 5.5 | 174 | 120.200.021 | FMS |
| Metal hub with gold surface | Reusable instruments into sterile packaging with 1 piece inside | | | |
| DR WAVE CUTTER | Ø 4.0 | 130 | 020.101.528 | FMS |
| DR WAVE CUTTER | Ø 5.0 | 130 | 020.201.528 | FMS |
| DR PIT BULL | Ø 4.0 | 130 | 020.101.538 | FMS |
| DR PIT BULL | Ø 5.0 | 130 | 020.201.538 | FMS |
| DR PIT BULL | Ø 7.0 | 130 | 020.701.538 | FMS |
| DR PIT BULL | Ø 4.0 | 155 | 020.101.535 | FMS |
| DR PIT BULL | Ø 5.0 | 155 | 020.201.535 | FMS |
| DR SYNOVIA PIT BULL | Ø 4.0 | 130 | 020.101.534 | FMS |
| DR BARREL BUR | Ø 4.0 | 130 | 130.100.519 | FMS |
| DR BARREL BUR | Ø 5.5 | 130 | 130.200.519 | FMS |
| DR ROUND BUR | Ø 4.0 | 130 | 120.100.519 | FMS |
| DR ROUND BUR | Ø 5.5 | 130 | 120.200.519 | FMS |

Gold: Can be used 10 times



DR WAVE CUTTER®
(one double blade)




DR PIT BULL®
(three double blades, large window)



DR SYNOVIA PIT BULL®
(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our FMS DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by JOHNSON & JOHNSON. The firm Johnson & Johnson Corporation Inc. and DR MEDICAL AG are independent companies.

Version001_EN