

DR Shaver Blades and Burs for the Johnson & Johnson FMS VUE Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	024.801.230	DRNX
DR WAVE CUTTER	Ø 3.0	85	024.301.231	DRNX
DR WAVE CUTTER	Ø 3.0	130	024.301.230	DRNX
DR WAVE CUTTER	Ø 4.0	130	024.101.230	DRNX
DR WAVE CUTTER	Ø 5.0	130	024.201.230	DRNX
DR PIT BULL	Ø 3.0	130	024.301.240	DRNX
DR PIT BULL	Ø 4.0	130	024.101.240	DRNX
DR PIT BULL	Ø 5.0	130	024.201.240	DRNX
DR PIT BULL	Ø 7.0	130	024.701.240	DRNX
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	024.101.236	DRNX
DR BENDING PIT BULL	Ø 3.0	130	024.301.340	DRNX
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	024.101.241	DRNX
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	024.101.242	DRNX
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	024.101.250	DRNX
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	024.101.251	DRNX
DR BARREL BUR	Ø 4.0	130	134.100.220	DRNX
DR BARREL BUR	Ø 5.5	130	134.200.220	DRNX
DR ROUND BUR	Ø 4.0	130	124.100.220	DRNX
DR ROUND BUR	Ø 5.5	130	124.200.220	DRNX
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	024.801.030	DRNX
DR WAVE CUTTER	Ø 3.0	85	024.301.030	DRNX
DR WAVE CUTTER	Ø 3.0	130	024.301.031	DRNX
DR WAVE CUTTER	Ø 4.0	130	024.101.030	DRNX
DR WAVE CUTTER	Ø 5.0	130	024.201.030	DRNX
DR PIT BULL	Ø 3.0	85	024.301.047	DRNX
DR PIT BULL	Ø 3.0	130	024.301.040	DRNX
DR PIT BULL	Ø 4.0	130	024.101.040	DRNX
DR PIT BULL	Ø 5.0	130	024.201.040	DRNX
DR PIT BULL	Ø 7.0	130	024.701.040	DRNX
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	024.101.036	DRNX
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	024.201.043	DRNX
DR PIT BULL	Ø 4.0	155	024.101.037	DRNX
DR PIT BULL	Ø 5.0	155	024.201.037	DRNX
DR BENDING PIT BULL	Ø 3.0	130	024.301.140	DRNX
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	024.101.041	DRNX
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	024.101.042	DRNX
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	024.101.050	DRNX
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	024.101.051	DRNX
DR BARREL BUR	Ø 4.0	130	134.100.020	DRNX
DR BARREL BUR	Ø 5.5	130	134.200.020	DRNX
DR BARREL BUR	Ø 4.0	174	134.100.021	DRNX
DR ROUND BUR	Ø 4.0	130	124.100.020	DRNX
DR ROUND BUR	Ø 5.5	130	124.200.020	DRNX
DR ROUND BUR	Ø 4.0	174	124.100.021	DRNX
DR ROUND BUR	Ø 5.5	174	124.200.021	DRNX
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	024.101.528	DRNX
DR WAVE CUTTER	Ø 5.0	130	024.201.528	DRNX
DR PIT BULL	Ø 4.0	130	024.101.538	DRNX
DR PIT BULL	Ø 5.0	130	024.201.538	DRNX
DR PIT BULL	Ø 7.0	130	024.701.538	DRNX
DR PIT BULL	Ø 4.0	155	024.101.535	DRNX
DR PIT BULL	Ø 5.0	155	024.201.535	DRNX
DR SYNOVIA PIT BULL	Ø 4.0	130	024.101.534	DRNX
DR BARREL BUR	Ø 4.0	130	134.100.519	DRNX
DR BARREL BUR	Ø 5.5	130	134.200.519	DRNX
DR ROUND BUR	Ø 4.0	130	124.100.519	DRNX
DR ROUND BUR	Ø 5.5	130	124.200.519	DRNX

Gold: Can be used 10 times



DR WAVE CUTTER®

(one double blade)



DR PIT BULL®

(three double blades, large window)



DR SYNOVIA PIT BULL®

(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRNX DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by JOHNSON & JOHNSON. The firm Johnson & Johnson Corporation Inc. and DR MEDICAL AG are independent companies.

Version001_EN