

DR Shaver Blades and Burs for the DYONICS Smith & Nephew Handpiece



Our DR SHAVER BLADES AND BURS distinguish themselves with:

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

DR WAVE CUTTER®

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

DR PIT BULL®

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

DR SHAVER BUR CERATITAN®

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made !

DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
Plastic blue hub	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	021.801.230	DRSN
DR WAVE CUTTER	Ø 3.0	85	021.301.231	DRSN
DR WAVE CUTTER	Ø 3.0	130	021.301.230	DRSN
DR WAVE CUTTER	Ø 4.0	130	021.101.230	DRSN
DR WAVE CUTTER	Ø 5.0	130	021.201.230	DRSN
DR PIT BULL	Ø 3.0	130	021.301.240	DRSN
DR PIT BULL	Ø 4.0	130	021.101.240	DRSN
DR PIT BULL	Ø 5.0	130	021.201.240	DRSN
DR PIT BULL	Ø 7.0	130	021.701.240	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.236	DRSN
DR BENDING PIT BULL	Ø 3.0	130	021.301.340	DRSN
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.241	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.242	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.250	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.251	DRSN
DR BARREL BUR	Ø 4.0	130	131.100.220	DRSN
DR BARREL BUR	Ø 5.5	130	131.200.220	DRSN
DR ROUND BUR	Ø 4.0	130	121.100.220	DRSN
DR ROUND BUR	Ø 5.5	130	121.200.220	DRSN
Metal hub silver	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	021.801.030	DRSN
DR WAVE CUTTER	Ø 3.0	85	021.301.030	DRSN
DR WAVE CUTTER	Ø 3.0	130	021.301.031	DRSN
DR WAVE CUTTER	Ø 4.0	130	021.101.030	DRSN
DR WAVE CUTTER	Ø 5.0	130	021.201.030	DRSN
DR PIT BULL	Ø 3.0	85	021.301.047	DRSN
DR PIT BULL	Ø 3.0	130	021.301.040	DRSN
DR PIT BULL	Ø 4.0	130	021.101.040	DRSN
DR PIT BULL	Ø 5.0	130	021.201.040	DRSN
DR PIT BULL	Ø 7.0	130	021.701.040	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.036	DRSN
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	021.201.043	DRSN
DR PIT BULL	Ø 4.0	155	021.101.037	DRSN
DR PIT BULL	Ø 5.0	155	021.201.037	DRSN
DR BENDING PIT BULL	Ø 3.0	130	021.301.140	DRSN
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.041	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	021.101.042	DRSN
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.050	DRSN
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	021.101.051	DRSN
DR BARREL BUR	Ø 4.0	130	131.100.020	DRSN
DR BARREL BUR	Ø 5.5	130	131.200.020	DRSN
DR BARREL BUR	Ø 4.0	174	131.100.021	DRSN
DR ROUND BUR	Ø 4.0	130	121.100.020	DRSN
DR ROUND BUR	Ø 5.5	130	121.200.020	DRSN
DR ROUND BUR	Ø 4.0	174	121.100.021	DRSN
DR ROUND BUR	Ø 5.5	174	121.200.021	DRSN
Metal hub with gold surface	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	021.101.528	DRSN
DR WAVE CUTTER	Ø 5.0	130	021.201.528	DRSN
DR PIT BULL	Ø 4.0	130	021.101.538	DRSN
DR PIT BULL	Ø 5.0	130	021.201.538	DRSN
DR PIT BULL	Ø 7.0	130	021.701.538	DRSN
DR PIT BULL	Ø 4.0	155	021.101.535	DRSN
DR PIT BULL	Ø 5.0	155	021.201.535	DRSN
DR SYNOVIA PIT BULL	Ø 4.0	130	021.101.534	DRSN
DR BARREL BUR	Ø 4.0	130	131.100.519	DRSN
DR BARREL BUR	Ø 5.5	130	131.200.519	DRSN
DR ROUND BUR	Ø 4.0	130	121.100.519	DRSN
DR ROUND BUR	Ø 5.5	130	121.200.519	DRSN

Gold: Can be used 10 times



DR WAVE CUTTER®
(one double blade)



DR PIT BULL®
(three double blades, large window)



DR SYNOVIA PIT BULL®
(three double blades, small window)



DR BARREL BUR®



DR ROUND BUR®

Our DRSN DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by SMITH & NEPHEW. The firm Smith & Nephew Corporation Inc. and DR MEDICAL AG are independent companies.

Version001_EN